




Endoscopic Accessories Product Catalog



PROFESSIONAL MANUFACTURER FOR
PREMIUM HEALTHCARE PRODUCTS



VeriEndo® Single Use Biopsy Forceps



*Reliable and Easy Cut
Get Ideal Samples*

Used to obtain tissue sample for the diagnosis of suspected pathology with an endoscope

- Premium medical-grade stainless material
- Sharp cup edge to get samples with less bleeding and quickly
- Rigid four-link structure makes the cup opened and closed smoothly and easily
- Smooth surface minimizes damage to the endoscope
- Various specifications are available for your choice

 Oval Cup Biopsy Forceps	REF	Working Length (mm)	Outer Diameter (mm)	Minimum Channel (mm)
	EBF-A-1806 ★	600	1.8	2.0
	EBF-A-1812 ★	1200	1.8	2.0
	EBF-A-1816	1600	1.8	2.0
	EBF-A-1818	1800	1.8	2.0
	EBF-A-2412	1200	2.4	2.8
	EBF-A-2416 ★	1600	2.4	2.8
	EBF-A-2418 ★	1800	2.4	2.8
	EBF-A-2423 ★	2300	2.4	2.8

 Oval Cup Biopsy Forceps with Spikes	REF	Working Length (mm)	Outer Diameter (mm)	Minimum Channel (mm)
	EBF-B-1812	1200	1.8	2.0
	EBF-B-1816	1600	1.8	2.0
	EBF-B-1818	1800	1.8	2.0
	EBF-B-2412	1200	2.4	2.8
	EBF-B-2416 ★	1600	2.4	2.8
	EBF-B-2418 ★	1800	2.4	2.8
	EBF-B-2423 ★	2300	2.4	2.8

 Alligator Cup Biopsy Forceps	REF	Working Length (mm)	Outer Diameter (mm)	Minimum Channel (mm)
	EBF-C-1112	1200	1.1	1.2
	EBF-C-1806 ★	600	1.8	2.0
	EBF-C-1812 ★	1200	1.8	2.0
	EBF-C-1816	1600	1.8	2.0
	EBF-C-1818	1800	1.8	2.0
	EBF-C-2412	1200	2.4	2.8
	EBF-C-2416 ★	1600	2.4	2.8
	EBF-C-2418 ★	1800	2.4	2.8
	EBF-C-2423 ★	2300	2.4	2.8

 Alligator Cup Biopsy Forceps with Spikes	REF	Working Length (mm)	Outer Diameter (mm)	Minimum Channel (mm)
	EBF-D-2412	1200	2.4	2.8
	EBF-D-2416 ★	1600	2.4	2.8
	EBF-D-2418 ★	1800	2.4	2.8
	EBF-D-2423 ★	2300	2.4	2.8

※ Various working lengths and cup diameters are available. For more specifications, please contact us at info@promisemed.ca

PROMISEMED HANGZHOU MEDITECH CO., LTD.


No. 1388 Cangxing Street, Cangqian Community,
Yuhang District, Hangzhou City, 311121 Zhejiang, China
info@promisemed.ca
www.promisemed-global.com



6666-EBF-VeriEndo-H-02 (EN)



VeriEndo® Single Use Cytology Brushes



Used to get tissue sample for the diagnosis of suspected pathology with an endoscope

- Thinner bristles evenly distributed in large spread angle for optimal sampling
- Ideal brush hardness reduces tissue trauma
- The brush head can be advanced and withdrawn smoothly
- Specially processed brush head reduces damage to the endoscope and the patient
- Quick exchange cytology brushes (Type DL) used for biliary tract brushing

 Type L	REF	Working Length (mm)	Max. Diameter of the Insertion Part (mm)	Brush Length (mm)	Brush Diameter (mm)
	ECB-0720-1020-L	700	2.0	10	2.0
	ECB-1220-1012-L	1200	2.0	10	1.2
	ECB-1220-1020-L ★	1200	2.0	10	2.0
	ECB-1220-1030-L ★	1200	2.0	10	3.0
	ECB-1620-1020-L	1600	2.0	10	2.0
	ECB-2020-1020-L	2000	2.0	10	2.0
	ECB-2320-1020-L	2300	2.0	10	2.0
	ECB-1628-1030-L ★	1600	2.8	10	3.0
	ECB-2028-1030-L	2000	2.8	10	3.0
	ECB-2328-1030-L ★	2300	2.8	10	3.0

 Type U	REF	Working Length (mm)	Max. Diameter of the Insertion Part (mm)	Brush Length (mm)	Brush Diameter (mm)
	ECB-1628-0860-U	1600	2.8	8	6.0
	ECB-2028-0860-U	2000	2.8	8	6.0
	ECB-2328-0860-U	2300	2.8	8	6.0

 Type DL	REF	Working Length (mm)	Max. Diameter of the Insertion Part (mm)	Brush Length (mm)	Brush Diameter (mm)
	ECB-2320-1020-DL	2300	2.0	10	2.0
	ECB-2328-1030-DL ★	2300	2.8	10	3.0

※ For more specifications, please contact us at info@promisemed.ca

PROMISEMED HANGZHOU MEDITECH CO., LTD.


No. 1388 Cangxing Street, Cangqian Community,
Yuhang District, Hangzhou City, 311121 Zhejiang, China
info@promisemed.ca
www.promisemed-global.com



6666-EBF-VeriEndo-H-02 (EN)



VeriEndo® Single Use Grasping Forceps



Used to take out foreign bodies or excised tissue with an endoscope

- Large opening width makes it easy to catch large foreign bodies
- The cup can be opened and closed smoothly and easily
- Good gripping performance
- Flexible and coated sheath to aid smooth insertion to an endoscope
- Various specifications are available for your choice

	REF	Working Length (mm)	Outer Diameter (mm)	Minimum Channel (mm)
Alligator Grasping Forceps 	EGF-A-2406-N	600	2.4	2.8
	EGF-A-2412-N	1200	2.4	2.8
	EGF-A-2418-N ★	1800	2.4	2.8
	EGF-A-2423-N ★	2300	2.4	2.8
	EGF-A-2427-N	2700	2.4	2.8

	REF	Working Length (mm)	Outer Diameter (mm)	Minimum Channel (mm)
Rat Tooth Grasping Forceps 	EGF-B-2406-N	600	2.4	2.8
	EGF-B-2412-N	1200	2.4	2.8
	EGF-B-2418-N ★	1800	2.4	2.8
	EGF-B-2423-N ★	2300	2.4	2.8
	EGF-B-2427-N	2700	2.4	2.8

	REF	Working Length (mm)	Outer Diameter (mm)	Minimum Channel (mm)
Rat Tooth-Alligator Grasping Forceps 	EGF-C-2406-N	600	2.4	2.8
	EGF-C-2412-N	1200	2.4	2.8
	EGF-C-2418-N	1800	2.4	2.8
	EGF-C-2423-N ★	2300	2.4	2.8
	EGF-C-2427-N	2700	2.4	2.8

※ For more specifications, please contact us at info@promisemed.ca

PROMISEMED HANGZHOU MEDITECH CO., LTD.


No. 1388 Cangxing Street, Cangqian Community,
Yuhang District, Hangzhou City, 311121 Zhejiang, China
info@promisemed.ca
www.promisemed-global.com



6666-EBF-VeriEndo-H-01 (EN)



VeriEndo® Retrieval Nets



Used to take out foreign bodies or excised tissue with an endoscope

- 360°rotation for flexible position
- Strong net to greatly reduce tearing
- The net can be opened and closed smoothly and easily
- PTFE sheath with the feature of high lubrication for smooth insertion

REF	Sheath O.D. (mm)	Working Length (mm)	Working Width (mm)
-----	---------------------	------------------------	-----------------------



ERN-H1815A ★	2.4	1800	15
ERN-H1825A ★	2.4	1800	25
ERN-H1832A	2.4	1800	32
ERN-H2315A ★	2.4	2300	15
ERN-H2325A ★	2.4	2300	25
ERN-H2332A	2.4	2300	32
ERN-H2715A	2.4	2700	15
ERN-H2725A	2.4	2700	25
ERN-H2732A	2.4	2700	32

PROMISEMED HANGZHOU MEDITECH CO., LTD.

No. 1388 Cangxing Street, Cangqian Community,
Yuhang District, Hangzhou City, 311121 Zhejiang, China
info@promisemed.ca
www.promisemed-global.com





VeriEndo[®] 3 Prongs Grasping Forceps



Used to take out foreign bodies or excised tissue with an endoscope

- Hooked prongs promote secure grip and reliable retrieval
- The prongs can be advanced and retracted smoothly and easily
- PTFE sheath with the feature of high lubrication for smooth insertion and is not easy to be deformed

REF	Sheath O.D. (mm)	Working Length (mm)
ERN-H1830B ★	2.4	1800
ERN-H2330B ★	2.4	2300
ERN-H2730B ★	2.4	2700

PROMISEMED HANGZHOU MEDITECH CO., LTD.


No. 1388 Cangxing Street, Cangqian Community,
Yuhang District, Hangzhou City, 311121 Zhejiang, China
info@promisemed.ca
www.promisemed-global.com



6666-ERN-VeriEndo-B-H-02 (EN)




VeriEndo[®] Endoscopic Injection Needles



Used to inject liquid into gastrointestinal mucosa with an endoscope

- Unique needle bevel allows small penetration force, ensuring good puncture
- Special shaped inner tube dramatically avoids blockage and coiling
- Optimized inner tube structural size reduces injection resistance and increases flow rate
- Safety lock design avoids the needle being accidentally advanced during insertion to damage the endoscopy
- Ergonomic handle for easy and comfortable operation

	REF	Sheath O.D. (mm)	Working Length (mm)	Needle Gauge	Needle Length (mm)
	EIN-H182104	2.4	1800	21G (0.8 mm)	4
	EIN-H182106	2.4	1800	21G (0.8 mm)	6
	EIN-H182304	2.4	1800	23G (0.6 mm)	4
	EIN-H182306	2.4	1800	23G (0.6 mm)	6
	EIN-H182504	2.4	1800	25G (0.5 mm)	4
	EIN-H182506	2.4	1800	25G (0.5 mm)	6
	EIN-H242104	2.4	2300	21G (0.8 mm)	4
	EIN-H242106	2.4	2300	21G (0.8 mm)	6
	EIN-H242304	2.4	2300	23G (0.6 mm)	4
	EIN-H242306	2.4	2300	23G (0.6 mm)	6
	EIN-H242504	2.4	2300	25G (0.5 mm)	4
	EIN-H242506	2.4	2300	25G (0.5 mm)	6

※ Various working lengths and cup diameters are available. For more specifications, please contact us at info@promisemed.ca

PROMISEMED HANGZHOU MEDITECH CO., LTD.

No. 1388 Cangxing Street, Cangqian Community,
Yuhang District, Hangzhou City, 311121 Zhejiang, China
info@promisemed.ca
www.promisemed-global.com



6666-EIN-VeriEndo-HI-01 (EN)



VeriEndo® Polypectomy Snares




Used to resect polyps (recommended size of 10-19mm) with an endoscope

- 360°rotation with stable rotation performance for flexible resection
- Special braided wire with good cutting and blood coagulation performance
- Premium sheath material with good deforming resistance
- The snare can be opened and closed easily and smoothly
- Various specifications for your choice



Oval

REF	Sheath O.D. (mm)	Working Length (mm)	Loop Width (mm)
ESN-S1810T	1.8	1800	10
ESN-S1815T	1.8	1800	15
ESN-H1810T	2.4	1800	10
ESN-H1815T ★	2.4	1800	15
ESN-H1825T	2.4	1800	25
ESN-H1832T	2.4	1800	32
ESN-H2410T	2.4	2400	10
ESN-H2415T ★	2.4	2400	15
ESN-H2425T	2.4	2400	25
ESN-H2432T	2.4	2400	32

	REF	Sheath O.D. (mm)	Working Length (mm)	Loop Width (mm)
	ESN-H2410L	2.4	2400	10
	ESN-H2415L	2.4	2400	15
	ESN-H2422L	2.4	2400	22
	ESN-H2432L	2.4	2400	32
	Hexagonal			

 Crescent	REF	Sheath O.D. (mm)	Working Length (mm)	Loop Width (mm)
	ESN-H2415B	2.4	2400	15
	ESN-H2425B	2.4	2400	25

※ For more specifications, please contact us at info@promisemed.ca

PROMISEMED HANGZHOU MEDITECH CO., LTD.

No. 1388 Cangxing Street, Cangqian Community,
Yuhang District, Hangzhou City, 311121 Zhejiang, China
info@promisemed.ca
www.promisemed-global.com



6666-ESN-VeriEndo-H-00 (EN)



VeriEndo® Cold Snares



Used to resect diminutive polyps (<5 mm) or noncancerous polyps up to 10 mm in size with an endoscope

- Thin wire to reduce the risk of specimen flying away from the resection site
- 360° rotation with stable rotation performance for flexible resection
- Special designed wire with good cutting performance ensures that polyps can be resected cleanly
- Stiffer sheath aids in reducing buckling during resection
- Ergonomic handle for easy and comfortable operation

REF	Sheath O.D. (mm)	Working Length (mm)	Loop Width (mm)
ECS-H1610T	2.4	1600	10
ECS-H1615T	2.4	1600	15
ECS-H1625T	2.4	1600	25
ECS-H1630T	2.4	1600	32
ECS-H1810T	2.4	1800	10
ECS-H1815T	2.4	1800	15
ECS-H1825T	2.4	1800	25
ECS-H1830T	2.4	1800	32
ECS-H2310T ★	2.4	2300	10
ECS-H2315T ★	2.4	2300	15
ECS-H2325T	2.4	2300	25
ECS-H2332T	2.4	2300	32
ECS-H2710T	2.4	2700	10
ECS-H2715T	2.4	2700	15
ECS-H2725T	2.4	2700	25
ECS-H2732T	2.4	2700	32

※ For more specifications, please contact us at info@promisemed.ca

PROMISEMED HANGZHOU MEDITECH CO., LTD.

No. 1388 Cangxing Street, Cangqian Community,
Yuhang District, Hangzhou City, 311121 Zhejiang, China
info@promisemed.ca
www.promisemed-global.com



6666-EBF-VeriEndo-H-02 (EN)



VeriEndo[®] Spraying Catheter for Single Use



Used to spray stain or solution with an endoscope

- Improved nozzle for fine mist, large spraying angle and even distribution
- Special shaped inner tube dramatically avoids blockage and coiling
- PTFE sheath with the feature of high lubrication and slippery spray tip for smooth insertion

PROMISEMED HANGZHOU MEDITECH CO., LTD.


No. 1388 Cangxing Street, Cangqian Community,
Yuhang District, Hangzhou City, 311121 Zhejiang, China
info@promisemed.ca
www.promisemed-global.com



6666-EIN-VeriEndo-H-00 (EN)



VeriEndo[®] Single Use Hemoclip



Used to be placed in the gastrointestinal tract with an endoscope for the purpose of hemostasis, endoscopic marking and so on

- Repeatability of opening and closing attribute to achieve repositioning
- 360° and 1:1 rotation for flexible position
- Strong clip retention
- Ergonomic handle for easy and comfortable operation
- Multiple opening widths meet various needs

	REF	Working Length (mm)	Opening Width (mm)	Claw Angle	Outer Diameter (mm)
	HRA-27012-135-Y	2700	12	135°	2.6
	HRA-23012-135-Y ★	2300	12	135°	2.6
	HRA-19512-135-Y ★	1950	12	135°	2.6
	HRA-16512-135-Y	1650	12	135°	2.6

※ For more specifications, please contact us at info@promisemed.ca

PROMISEMED HANGZHOU MEDITECH CO., LTD.

No. 1388 Cangxing Street, Cangqian Community,
Yuhang District, Hangzhou City, 311121 Zhejiang, China
info@promisemed.ca
www.promisemed-global.com

